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Envirobeads S-X3

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GPC

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PCB

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| | Poultry Lipid | Pork Lipid | Beef Lipid |
|---------------------|---------------|------------|------------|
| Average %Recovery | 98.3 | 96.5 | 98.0 |
| %Standard Deviation | 3.7 | 5.0 | 6.0 |

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1. Ault, J.A. ; Schofield, M.S. ; Johnson, L.D. ; Waltz, R.H. Automated Gel Permeation Chromatographic Preparation of Vegetables, Fruits, and Crops for Organophosphate Residues Determination Utilizing Flame Photometric Detection. *J. Agricultural and Food Chem.* 1979, *27*, 825.
2. Fehringer, N.V. Determination of Polybrominated Biphenyl Residues in Dairy Products. *J. A. O. A. C.* 1975, *58(5)*, 978.
3. Hopper, M.L. Evaluation of Gel Permeation Chromatography (GPC) as a General Cleanup Tool in FDA' s Multiresidue Methodology. FDA Kansas City District, Final Report, SARP No. 78 - 76, 1979.
4. Johnson, L.D. ; Waltz, R.H. ; Ussary, J.P. ; Kaiser, F.E. Automated Gel Permeation Chromatographic Cleanup of Animal and Plant Extracts for Pesticide Residue Determination. *J. A. O. A. C.* 1976, *59(1)*, 174.
5. Moseman, R.F. ; Ward, M.K. ; Christ, H.L. ; Zehr, R.D. Kepone in Human Blood and Environmental Samples. *J. Agri. Food Chem.* 1978, *26*, 965.
6. Reickel, W.L. ; Kolbe, E. ; Stafford, C.J. Gas-Liquid Chromatographic and Gas-Liquid Chromatographic-Mass Spectrometric Determination of Fenvalerate and Permethrin Residues in Grasshoppers and Duck Tissue Samples. *J. A. O. A. C.* 1981, *64(5)*, 1196 - 1200.
7. Stalling, D.L. ; Johnson, J. ; Huckens, J.N. Automated Gel Permeation-Carbon Chromatographic Cleanup of Dioxins, PCB' s, Pesticides, and Industrial Chemicals. *Environ. Qual. and Safety Supplement, Volume III*, IUPAC Third International Congress of Pesticide Chemistry, Helsinki, 3 - 9 July, 1974.
8. Stalling, D.L. ; Tindle, R.C. ; Johnson, J.L. Cleanup of Pesticide and Polychlorinated Biphenyl Residues in Fish Extracts by Gel Permeation Chromatography. *J. A. O. A. C.* 1972, *55(1)*, 32.

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